



CORNING
Discovering Beyond Imagination



**Life
Sciences**

SnAPPShots

*A brief report
from the Corning
Applications Group*

Processing Arrays with Pronto!™ Universal Hybridization System in Tecan HS 4800 Hybridization Station


Camilo Canel, Michael W. Briggs, Laurent Picard, and Debra S. Hoover

The reliability of microarray data is affected by a large number of factors. Foremost is the quality of the substrate and the reagents used to manufacture and process the microarrays.

Further improvement in reproducibility and sensitivity may be achieved through automation. Researchers have long sought the ability to perform multiple array hybridizations under tightly controlled and reproducible conditions. Array experiments typically call for processing at least four arrays in parallel, which, when carried out in single-slide hybridization chambers and glass containers, requires significant operator intervention and introduces a high level of experimental error that decreases the sensitivity of the assays.


Corning welcomes the availability of reliable automated systems and seeks to facilitate their use by providing protocols for the use of our market-leading Pronto! hybridization reagents in these automated environments.





Highly Reproducible Hybridization of cDNAs Arrayed on Corning® UltraGAPS™ Slides

Four human cDNA arrays obtained from the Institute for Genomic Research (TIGR) were treated with the background-reduction reagents in the Pronto! Background Reduction Kit (Cat. No. 40029) in glass containers, dried, and placed in a Tecan HS 4800 hybridization station. The other steps, from pre-hybridization to drying, were completed using Corning Pronto! hybridization reagents exclusively and Pronto! Plus labeling reagents exclusively. The CV's across the arrays tested were 15% in both channels.



The protocol outlined above has been optimized for hybridization of fluorescently labeled cDNA (transcriptional profiling) and genomic DNA (comparative genomic hybridization) to arrays of cDNA and long oligonucleotides, using the Corning® Pronto!™ Universal Hybridization Kit. The protocol includes steps for both prehybridization and hybridization.

For further information about the Pronto! reagents, please visit www.corning.com/lifesciences. Contact Dr. Camilo Canel (207.985.5353, CanelC@corning.com), or, in Europe, Dr. Laurent Picard (31.20.655.7942, PicardL@corning.com), to request a program file (.hpr) to run this protocol on the Tecan hybridization stations. For information about the Tecan HS 4800 System, visit www.tecan.com.

Corning Incorporated Life Sciences

45 Nagog Park
Acton, MA 01720
t 800.492.1110
t 978.635.2200
f 978.635.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA

Australia
t 61 2-9416-0492
f 61 2-9416-0493

China
t 86 21-3222-4666
f 86 21-6288-1575

Hong Kong
t 852-2807-2723
f 852-2807-2152

India

t 91 11 341 3440
f 91 11 341 1520

Japan

t 81 (0) 3-3586 1996/1997
f 81 (0) 3-3586 1291/1292

Korea

t 82 2-796-9500
f 82 2-796-9300

Singapore

t 65 6733-6511
f 65 6861-2913

Taiwan

t 886 2-2716-0338
f 886 2-2716-0339

EUROPE

France

t 0800 916 882
f 0800 918 636

Germany

t 0800 101 1153
f 0800 101 2427

The Netherlands

& All Other European Countries

t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

United Kingdom

t 0800 376 8660
f 0800 279 1117

LATIN AMERICA

Brasil

t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico

t (52-81) 8158-8400
f (52-81) 8313-8589