

## Soil respiration measurement using Micro-Oxymax Respirometer

Nov. 2008  
Sanghoon Kang

### Basal respiration

- I. Soil preparation
  - a. Pre-incubate several hours to warm up at incubation temperature (RT)<sup>1</sup>
  - b. Weight wet soil 50g dry weight based on moisture content measurement into numbered (matching to channel number) Duran bottles.
  - c. After all bottles are weighed, finish assembling with water bottle and sample bottles making sure O-rings and metal plates connects are intact.
- II. Setting up the Micro Oxymax Respirometer
  - a. Refill drierite in drying columns
    - i. Dry used drierite (pink) in drying oven, take dried drierite from the oven if there is drierite in the oven (color should be little bit bluish or at least purple).
  - b. Calibration (Tools-Calibration)
    - i. Offsetting gases other than O<sub>2</sub> (bottle 1 – lime soda column)
    - ii. Offsetting O<sub>2</sub> (bottle 2 – calibration gas cylinder<sup>2</sup>)
    - iii. Gaining O<sub>2</sub> (bottle 1)
    - iv. Gaining gases other than O<sub>2</sub> (bottle 2<sup>3</sup>)
  - c. Setup the experiment (Experiment-Setup)
    - i. Setup tab
      1. **Channels** – End channel – # channels to use
      2. Mark ‘O2 Consumption Positive’
      3. **Timing** – Sample Interval(Hours) – 2.5hours or ‘auto’
      4. **Data Units** – Normalization Unit – g
    - ii. Chamber Setup tab
      1. Type ‘Normalization Units’ and ‘Channel Label’ while leakage test
      2. ‘Volume’ and ‘Leakage’ will be determined by tests
  - d. Leakage check & volume measurement (Tools-Utilities)
    - i. Click on Leakage button for testing (pp. 76 for source of leakage)
    - ii. Click on Volume button for measuring headspace volume
  - e. Start the incubation by clicking Run

---

<sup>1</sup> Pre-incubation is needed for restoration of metabolic equilibrium of the population after storing at 4C (12hr) and re-wetting (24hr).

<sup>2</sup> Open the cylinder first, open regulator knob until hearing hissing sound, then adjust pressure to 5psi.

<sup>3</sup> Current composition (H<sub>2</sub>S 179.8ppm, H<sub>2</sub> 1801ppm, CO<sub>2</sub> 0.906%, CO 0.9015%, CH<sub>4</sub> 0.9054%) – balanced by N<sub>2</sub> (cylinder # SG9147282, reference # 83-124147696-1). It should be modified in System Properties at Sensors tab (pp. 70)

- III. Incubation condition & measurement
  - a. 23C (or RT), 2 ½ hour interval<sup>4</sup> for 120 hours<sup>5</sup>
  - b. Humidifying bottles (50ml with GL32 cap) for each sample to maintain moisture level in gas after going through drier (Figure 1)
  - c. Collect final results (accumulation) and graph in both rate (  $\mu\text{l O}_2 \text{ g}^{-1} \text{ soil h}^{-1}$ ) and accumulation ( $\mu\text{l O}_2 \text{ g}^{-1} \text{ soil}$ ) for further inference
- IV. After incubation
  - a. Finish incubation by clicking `Stop`.
  - b. Disconnect Duran bottles, empty soil, wash, oven-dry and cap with aluminum foil for long term storage; otherwise, re-use them for next incubation.
  - c. Leave water bottles connected for immediate next incubation; otherwise cover with the original orange caps

---

<sup>4</sup> Or shortest intervals possible. (# channel + 1) \* 6 min

<sup>5</sup> 72 hours by Chuck, 150 hours by Hollender *et al* 2003

# Experiment-Setup

**Experiment Running**

Setup Chamber Setup Comments Data Graphing

**Channels**  
Start Chamber: 1  
End Chamber: 6

**Timing**  
Sample Interval(Hours): 2.50 (2:150)  
Experiment Duration(Hours): N.A.

**Refresh**  
Refresh Threshold(%): 0.50  
Refresh interval: N.A.  
Refresh Window(sec): 38

**Data Units**  
Gas Units: ul  
Time Units: min  
Normalization Units: g

**Misc Setup**  
Auto Volume Measurement: ☒  
Purge Sensors: ☒  
Switch Driers: ☒  
O2 Consumption Positive: ☒  
Enable Open Flow Mode: ☐  
Manually Enter Chamber Temps.: ☐  
Aux. Temp Start Channel: N.A.

Open Stop Run Save Print Exit

Filename: \\data\\real test\\080918.dat Status: Experiment Running! Sys Waiting 00:48:18

**Experiment Running**

Setup **Chamber Setup** Comments Data Graphing

Ch	Volume (ml)	Normalize Units	Restriction (mmHg)	Leakage (ml/min)	Ch Temp (°C)	Channel Label
Ch 1	326	1,000 g		-0.080	Auto	blank
Ch 2	306	1,000 g		-0.058	Auto	positive control
Ch 3	202	20,000 g		0.064	Auto	field
Ch 4	198	20,000 g		0.114	Auto	50%
Ch 5	370	50,000 g		-0.010	Auto	field
Ch 6	350	50,000 g		-0.013	Auto	50%

Utilities

Open Stop Run Save Print Exit

Filename: \\data\\real test\\080918.dat Status: Experiment Running! Sys Waiting 00:49:10

